



Reliable

MEDRAD[®] Spectris Solaris EP Injection System

Reliable, efficient, patient focused.

Reliable

- Cleared for use with MR Scanners up to 3.0T
- Mobile configuration available for MR trailers
- Intuitive user controls
- Protocol memory saves time



Click & hold to rotate



Workflow efficiency

- Retrieve stored protocols
- Quick and easy to program
- Intuitive touchscreen
- Fast and easy to fill and prime
- Easy to install and remove syringes
- Flexible patient access



Patient Care

- Clear plastic syringes and tubing
- Easy to see possible air bubbles
- KVO delivers small boluses of fluid reducing the possibility of IV line clotting
- KVO runs independent from programmed injection

- Programmable pressure limit
- Real-time pressure display during injection